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## Literatur

- 1 Scholz J. (2008): Lungenprotektive Beatmung – Beatmungsinduzierte Lungenschädigung minimieren. *Anästhesiologie Intensivmedizin Notfallmedizin Schmerztherapie* 2008; 6: 436–437
- 2 Dreyfuss D, Soler P, Basset G, Saumon G (1988): High inflation pressure pulmonary edema. Respective effects of high airway pressure, high tidal volume and positive end-expiratory pressure. *American Review of Respiratory Disease* 1988; 137: 1159-1164
- 3 Kolobow T. (1988): Acute respiratory failure: On how to injure healthy lungs (and prevent sick lungs from recovering). *American Society for Artificial Internal Organs Transactions*. 1988; 34: 31-34
- 4 Kolobow T., Moretti MP, Fumagalli R, et al (1987): Serve impairment in the lung function included by high peak airway pressure during mechanical ventilation. *American Review of Respiratory Disease* 1988; 135: 312-315
- 5 The Acute Respiratory Distress Syndrome Network (2000): Ventilation with lower tidal volumes as compared with traditional tidal volumes for acute lung injury and the acute respiratory distress syndrome. *New England Journal of Medicine* 2000; 342: 130-1308
- 6 Dreyfuss D, Soler P, Basset G, Saumon G, (1988): High inflation pressure pulmonary edema. Respective effects of high airway pressure, high tidal volume, and positive end-expiratory pressure. *American Review of Respiratory Disease* 1988; 137: 1159–1164
- 7 Amato MB, Barbas CS, Medeiros DM et al (1998): Effect of a protective-ventilation strategy on mortality in the acute respiratory distress syndrome. *New England Journal of Medicine* 1998; 338:347–354
- 8 The Acute Respiratory Distress Syndrome Network (2000): Ventilation with lower tidal volumes as compared with traditional tidal volumes for acute lung injury and the acute respiratory distress syndrome. *New England Journal of Medicine* 2000; 342: 1301-1308

[Academy for respiratory therapy e.V.](http://www.academy-respiratory-therapy.com)  
[info@academy-respiratory-therapy.com](mailto:info@academy-respiratory-therapy.com)